MADByTE Early Adopter Installation Guide

Step 1: Install Anaconda

- Go to anaconda.org/downloads
- Select Anaconda for Windows Python 3 version
- Download the .exe installer
- Run the installer and use all default options
 - While it is installing, you will see an option for "Add Anaconda as my default Python 3.7" checked by default
 - Do that.
- When finished, run the program called "Anaconda Prompt"

Step 2: Install MADByTE

- Unzip the file called 'madbyte_master.zip' into your downloads folder.
 - o You will need to keep this folder here.
- Simple click the "madbyte_install.bat" script
 - o This will setup a local python environment with all the required packages
- You can now launch MADByTE with the "MADByTE.bat" script
- OPTIONALLY: you can copy a shortcut of the MADByTE.bat script to your desktop